



ERP Softwares Inc.

For All Web & Software Solutions.

Contact : Mr. Shirish Patel Mobile : +91 93212 87688

Email : info@erpsoftwares.com

Website : www.erpsoftwares.com

Quality Control Analyser

Statistical **C**ontrol Software v1.0
Quality

Tablets

Capsules - Cap & Body

Injection

Powder

Chemicals....

Compatible with All Balances

Quality Control Analyser is a statistical analysis software useful Capsule Manufacturing Companies, Tablet Manufacturing Companies, Chemical Laboratories, and many other industries.



Some of the companies where the customised software has been designed is mentioned below.

- 1) ACG World Wide. - Associated Capsules Ltd. (Kandivali - Mumbai)
- 2) INTAS Pharma (Ahemdabad)
- 3) IPCA - Laboratories Ltd (Silvasa)
- 4) Novarties





It is capable of acquiring measurement readings directly from the weighing balance and feeding the product into the weighing balance automatically through a auto feeder machine if any.

- This Software is compatible with all weighing balances and it is used to generate Data & Graphical Analysis of Items Representations. eg. Capsules - (Body / Cap).
- This is the advanced and sophisticated software designed for maximum performance to reduce the manpower and manual errors at the time of production.
- Reading can be taken as 50, 150, 210, 270.
- You can view whether the Weight Analysis falls within the Tolerance specifications.
- You can view the average Weight Analysis for every production batch. Shows you clearly the Upper Limit, Target Value & the Lower Limit.
- You can also view the Summary of weight with variability within a dip.
- Unlimited Data Record Storage Capacity. (Statistical Data can be available at any time in future)
- Data Export Facility.
- This is window based software designed for better accuracy higher performance less manpower error free operations.
- Administrative controls are provided that help you configure users of the package, their access restrictions, database back-up etc.

Statistical Quality Control Software. - [Weight Readings]

Batch Entry Form Details of Batches Automatic Readings Graphical Analysis Data Analysis

ACM No. 3

Batch No. BATCH004

Bath No. BATH0045

Cap / Body. CAP BODY

Operator Name. PANKAG

Size. OEL

Colour. GREEN

Target Value. 13 mg.

Tolerance Val. 10 % 7 %

Readings: 9 Bar 270

Date. 50 Capsules 04

Time. 9 Bar 7:02:23 PM

Weighing Details

Time Spend: Upper Limit > 13.91

Previous Weight: Target = 13

Current Weight: 0 Lower Limit < 11.7

Com Port Open: Tare Weight Actual Weight

Change Image [6735 X 4095]

Export to Excel Export to Word

WEIGH SAVE PRINT SEARCH UPDATE REFRESH ALL DELETE CAL EXIT

